



# OPTIMIZE TECHNOLOGIES

## Press Release

### OPTI-SOLV® Extreme High Pressure Pre-Column Hand Tight Filter System

Oregon City, OR – June 28, 2007 – Optimize Technologies, a leading provider of parts and consumables for HPLC is proud to introduce the OPTI-SOLV® Extreme High Pressure Pre-Column Hand Tight Filter System for extreme high pressure applications. The combination of low dispersion, high pressure capability and ease of use in the OPTI-SOLV® Extreme High Pressure Pre-Column Filter is an unbeatable choice for your column protection.

**OPTI-SOLV® Extreme High Pressure Pre-Column Filter** is an ideal choice for protecting columns with small particle size material employing ultra high pressure techniques. Such techniques analyze complex samples in the most demanding applications, which can cut short the life of these expensive columns.

**OPTI-SOLV® Extreme High Pressure Pre-Column Hand Tight Filters** help extend that life and protect your column investment, without sacrificing performance.



### **Advantages of the OPTI-SOLV® Extreme High Pressure Pre-Column Filter System**

- Tested and approved for use at 20,000 PSI
- Patented Auto Adjusting ZDV Column Connection
- Hand Tight filter replacement - NO TOOLS!
- Low Volume, low dispersion cartridges

15-04-03086	<b>OPTI-SOLV® EXP™</b>	Pre-Column Filter Holder	
15-04-03093	0.5 µm, <b>OPTI-SOLV® EXP™</b>	Replacement Filter	(5 pack)
15-04-03094	0.5 µm, <b>OPTI-SOLV® EXP™</b>	Replacement Filter	(10 pack)
15-04-03096	0.2 µm, <b>OPTI-SOLV® EXP™</b>	Replacement Filter	(5 pack)
15-04-03097	0.2 µm, <b>OPTI-SOLV® EXP™</b>	Replacement Filter	(10 pack)

### Your local Distributor:

Sequant GmbH, Heckenweg 6b, 45721 Haltern am See, Germany

Tel: +49(0)2364 50446-0

E-mail: [info@sequant.de](mailto:info@sequant.de)

Fax: +49(0)2364 50446-11

Website: [www.sequant.de](http://www.sequant.de)