

Save up to  
**40%**

## SiliaSep™ BT Flash Cartridges 75iS, 75iM & 75iL

### Specifications

<b>Products Number:</b>	FLH-R10030B-75iS (200g) FLH-R10030B-75iM (400g) FLH-R10030B-75iL (800g)
<b>Silica Characteristics:</b>	40-63 µm, 60 Å pH 6.5 - 7.5
<b>Maximum pressure:</b>	100 psi (7 bar)
<b>Flow rate range:</b>	100 - 250 mL/min
<b>Sample loading:</b>	0.1% to 10% w/w silica gel



### Description of SiliaSep BT Flash Cartridges

SiliaSep BT flash cartridges are designed to be used on Biotage™ FLASH 75™ system. These cartridges are manufactured under ISO quality standards to ensure consistent and reliable performance. Packed with the **UltraPure SiliaFlash** silica gel, they offer the following:

- High purity silica gel
- Neutral pH (6.5 - 7.5)
- Consistent water level content
- Innovative packing technology
- Leak-free guaranteed
- Lot-to-lot reproducibility
- High resolution
- No tailing

### SiliaSep BT Flash Cartridges (2/box)

Product Number	Silica
FLH-R10030B-75iS	200g
FLH-R10030B-75iM	400g
FLH-R10030B-75iL	800g

**Contact us for more information**